

## SPECIFICATION SHEET

**Product Name:** Human IgM **Catalog #:** L11010

**Description:** Human ImmunoglobulinM (IgM)

**Isotype:** IgM

**Source:** Normal Human Serum.

**Immunogen:** Purified Human IgM

**Concentration:** 2mg/ml **Lot:** L190401

**Buffer:** 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4

**Purification:** Octoic acid purified.

**Application:** Sepecific methodologies have not been tested using this product.

**Preservative:** None

**Storage& Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

**FOR RESEARCH USE ONLY**