

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #: L2G01201 **Lot #: L160601**

Description: Human Myeloperoxidase purified (MPO)

Grade: Affinity purified

Purity: 90% by SDS-PAGE

Source: Human leukocyte mass

Concentration: 3.3mg/ml

Form: 20mM Tris-HCl, pH 7.5, 0.5M NaCl

Biohazard: Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

Storage: -20°C or lower, for long term storage.