

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

### SPECIFICATION SHEET

<b>ProductName:</b>	OV-Ag	<b>Catalog#:</b>	L1F003
<b>Description:</b>	Human ovarian tissue Antigen ( OV-Ag )		
<b>Source:</b>	Human ovarian tissue		
<b>Format:</b>	Ion exchange, chromatography , Purified,Purified, MW: 12-120Kd Solubility is limited at neutral pH and physiological salt concentration.		
<b>Concentration&amp; Lot.:</b>	3.2mg/ml (Commassie Blue Staining Kit) <b>Lot:</b> L190301		
<b>Buffer:</b>	20mM Tris-Hcl 150Mm NaCl ,PH7.4		
<b>Preservative:</b>	0.2% NaN3		
<b>Applications:</b>	Suitable for use as a standard in immunoassay, Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp;Shelflife:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	A blood sample form the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		