

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Product Name:** IGFBP3 **Catalog #:** L4I00602

**Description:** Insulin Like Growth Factor Binding Protein 3(IGFBP3)

**Intorduction:**

Insulin-like growth factor-binding protein 3, also known as IGFBP-3, is a protein that in humans is encoded by the IGFBP3 gene. IGFBP-3 is one of six IGF binding proteins (IGFBP-1 to IGFBP-6) that have highly conserved structures and bind the insulin-like growth factors IGF-1 and IGF-2 with high affinity. IGFBP-7, sometimes inappropriately included in this family, shares neither the conserved structural features nor the high IGF affinity.

**Molecular weight:** 32KD

**Source:** E. coli

**Format:** Purified, Liquid

**Purification:** > 90% ,Sephrose derived purification

**Concentration& Lot.:** 1 mg/ml **Lot.:** L160701

**Buffer:** 20mm tris-hcl,300mM NaCl,0.2% SKL,PH 9.0

**Preservative:** 0.2% SKL

**Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw