

SPECIFICATION SHEET

Product Name: Re combinant cTnl protein **Catalog #:** L4C00102

Description: Human Troponin I-Cardiac >90% Recombinant, about 30KD

Source: Expressed in E. coli

Format: Purified

Purification: Purified, >90% pure by SDS-PAGE

Concentration&Lot: 1.1 mg/ml (Roche ECLIA) **Lot :** L20200901

Buffer: 10mmPBS, PH 7.4, 0.01%SKL

Preservative: None

Applications: Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.