

## SPECIFICATION SHEET

- Product Name:** Recombinant MYO Antigen      **Catalog #:** L4C00402
- Description:** Recombinant Human Myoglobin Antigen
- Source:** Expressed in E. coli
- Purification:** Purified, >85% pure by SDS-PAGE
- Concentration&Lot:** 1.5mg/ml (activity Concentration)      **Lot:** L20201201
- Buffer:** 50mM Tris-Hcl, 300mM NaCl, 0.1%SKL ,ph 8.0
- Preservative:** 0.1% SKL
- Applications:** Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
- Storage&Shelf life:** Store at 2~8°C for 3 years.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.