

## SPECIFICATION SHEET

<b>Product Name:</b>	Purified D-dimmer Antigen	<b>Catalog #:</b>	L4C00801
<b>Description:</b>	Purified native Human D-Dimer protein		
<b>Source:</b>	Human serum		
<b>Purification:</b>	Affinity purified.		
<b>Concentration&amp;Lot:</b>	5mg/ml	<b>Lot:</b>	L20200501
<b>Buffer:</b>	pH 7.5 PBS 150Mm NaCl		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.