

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	ICA 69	Catalog #:	L4I00403
Description:	Recombinant human Islet cells Antigen(ICA 69)		
Source:	Pichia pastoris , MW: about 72KD		
Format:	Purified, Liquid		
Purification:	> 85% (determined by SDS gelelectrophoresis)		
Concentration& Lot.:	1.66mg/ml	Lot.:	L190401
Buffer:	10mM TrisHCl,0.1%SKL,2mmEDTA		
Preservative:	0.1%SKL		
Application:	Immunoassays Immunization.		
Storage&Shelflife:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		