



Hainan Medikamende Trading Co., Ltd.

海南省迈迪卡门德商贸有限责任公司
Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	McAb to MCM5	Catalog #:	L1C03302				
Description:	McAb to Human MCM5 (Human Minichromosome Maintenance Deficient 5)Monoclonal Antibody						
Specificity:	Human Minichromosome Maintenance Deficient 5. No cross-reactivity with Human Albumin, Other cytokines.						
Clone:	MC002						
Host Animal:	Mouse. Hybridization of P3X63-Ag8.653 .myeloma cells with spleen cells from Balb/c mice.						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Purified MCM5 from Cord Blood						
Format:	Purified, Liquid						
Purification:	DEAE-chromatography						
Concentration&Lot:	Lot #:						
Buffer:	0.01M PB, 0.15M NaCl, pH 7.4.						
Preservative:	None						
Applications:	Determination of MCM5. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended antibody pair for immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1C03301</td><td>L1C03302</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1C03301	L1C03302
<u>Capture</u>	<u>Detection</u>						
L1C03301	L1C03302						
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.						

海南省迈迪卡门德商贸有限责任公司
地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
邮编：572000
电话：+86-898-31320087，传真：+86-898-31320087
网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.
Address: 6-A-502, Baoyuanhuayuan, Xiuyingdao-No.18,
Xiuyingqu, Haikou, 572000, Hainan, P.R. China
Phone: +86-898-31320087 Fax: +86-898-31320087
Website: www.medical-ivd.com Emai: info@medical-ivd.com