

Medical Electronic Detector New Product Digital Thermometer





Introduction

A thermometer is a device that measures temperature or a temperature gradient, has two important elements:

- (1) a temperature sensor in which some change occurs with a change in temperature;
- (2) some means of converting this change into a numerical value.

Thermometers are widely used in industry to monitor processes, in meteorology, in medicine, and in scientific research.



Quick Details

Type:	Thermometer, Rigid Tip
brand name	OEM
place of origin	Zhejiang, China
Model Number	MT-405
Instrument classification	class II
Product name	Medical Electronic Detector New Product Digital Thermometer
Certificate:	ISO 9001,ISO 13485,CE0197,RoHS,FDA
Key Words:	Pacificer,New Product Digital Thermometer Pacificer
MOQ:	2000 pcs
Battery type:	1.5V(LR41,SR41 or UCC 392)
LCD Display:	15.5*7mm
Unit Dimension:	6.2 cm x 4.2 cm x 4.5 cm (L x W x H)
Unit Weight:	Approx. 10 grams including battery
Warranty:	1year
SupplyAbility	3000000 Piece/Pieces per Month Medical Electronic Detector New Product Digital Thermometer



Product Description



Response Time	1 minute
Range	32.0°C – 42.9°C (90.0 °F – 109.9)
Accuracy	±0.1°C, 35.5°C – 42.0°C
	(±0.2°F, 95.9 °F-107.6 °F)
	±0.2°C below 35.5°C or over 42.0°C
	(±0.4 °F below 95.9 °F or over 107.6)
Display	Liquid crystal display, 4 1/2 digits,
Last Reading Memory, Fever Alarm, Auto-Shutoff	
Battery	One 1.5 V button battery included
	Size: LR41, SR41 or UCC 392
Battery life	About 2 years for average using



Detail Images





Related Products





Medical Electronic Detector New Product Digital Thermometer

New

Your Reliable Home Healthcare Medical Device Provider



- Digital Thermometer Specification:**
- Fast read
 - Waterproof
 - Backlight
 - Flexible tip
 - Jumbo LCD



- Blood Glucose System Specification:**
- Coding / no coding optional
 - Deluxe carry case
 - 5 seconds measuring time
 - 360 memories



Microprocessor single syringe infusion pump:

- self test
- Drop detection-optional
- Small size
- PC port-optional
- Large LCD
- 100 reading memories



Breast Pump Specification:

- Fits easily into purse or tote bag
- 2-phase expression technology
- All parts that contact with breast milk are made without BPA



- Hemoglobin System Specification:**
- Large LCD
 - 15 seconds measuring time
 - Deluxe carry case



- Blood Pressure Monitor Specification:**
- Jumbo LCD
 - 2+60 memories
 - Talking optional
 - Backlight optional
 - Irregular heartbeat detection

