

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	McAb to FFN	Catalog #:	L1F00101
Description:	Monoclonal Antibody to Fetal Fibronectin protein, FFN		
Clone:	FN001		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Fetal Fibronectin protein, FFN		
Format:	Purified, Liquid		
Concentration&Lot:	6.8 mg/ml (OD280nm) Lot #: L180601		
Buffer:	10mM PBS pH 7.4 ± 0.2		
Preservative:	None		
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay:		
Storage:	Store at -20°C for at least 3 years.		