

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

| | | | |
|-----------------------|---|---------------|---------|
| Catalog #: | L1V00101 | Lot #: | L190901 |
| Description: | Mouse anti-25-OH-Vitamin D3-KLH (Calcifediol-KLH) Monoclonal Antibody | | |
| Specificity: | Recognizes an epitope on the 25-OH-Vitamin D3-KLH(Calcifediol-KLH) | | |
| Immunogen: | 25-OH-Vitamin D3-KLH (Calcifediol-KLH) | | |
| Host Animal: | Mouse | | |
| Source: | Serum | | |
| Format: | Purified, Liquid | | |
| Concentration: | 1.8mg/ml | | |
| Buffer: | 0.01M PBS, pH 7.4 | | |
| Preservative: | None | | |
| Applications: | ELISA | | |
| Storage: | Store at -20°C for long term storage. Avoid repeated freeze-thaw cycles. | | |