



SPECIFICATION SHEET

Product Name:	McAb to Glycated hemoglobin	Catalog #:	L1K00504
Description:	Mouse anti-Glycated hemoglobin (hemoglobin A1c, HbA1c or HGBA1C) antibody		
Specificity:	Human		
Clone:	HB001		
Host Animal:	Mouse.		
Isotype:	IgG ₁		
Immunogen:	Purified Glycated hemoglobin (hemoglobin A1c, HbA1c or HGBA1C)		
Format:	Purified, Liquid		
Purification:	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4. No preservatives added.		
Concentration& Lot.:	4.6mg/ml	Lot. :	L190301
Buffer:	PBS, pH 7.4		
Applications:	Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		