

SPECIFICATION SHEET

| | | | |
|---------------------------------|--|-------------------|----------|
| Product Name: | Mouse anti human IgM-u chains | Catalog #: | L1100402 |
| Description: | Mouse Antibody to human Immunoglobulin M (IgM), u chains | | |
| Specificity: | Reacts with u chains on human IgM human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the u region against Normal human Serum. | | |
| Host Animal: | Mouse | | |
| Immunogen: | Purified Human IgM(u chains) | | |
| Format: | Ion exchange purification, Liquid | | |
| Purification: | >90% pure (SDS-PAGE). | | |
| Concentration: | 1.3 mg/ml | Lot.: | L170201 |
| Buffer: | 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4 | | |
| Preservative: | None | | |
| Storage& Shelf life: | Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles. | | |

FOR RESEARCH USE ONLY