

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Product Name:</b>	McAb to Human Albumin	<b>Catalog #:</b>	L3I00403
<b>Description:</b>	Mouse anti- Human Serum Albumin (HSA) antibody		
<b>Background:</b>	Human serum albumin is the most abundant protein in human blood plasma. It is produced in the liver. Albumin comprises about half of the blood serum protein. It is soluble and monomeric.		
<b>Specificity:</b>	Human		
<b>Clone:</b>	HA003		
<b>Host Animal:</b>	Mouse.		
<b>Isotype:</b>	IgG		
<b>Immunogen:</b>	Purified Human Serum Albumin		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4. No preservatives added.		
<b>Concentration&amp;Lot.:</b>	3.3 mg/ml (Lot. specific)		
<b>Buffer:</b>	PBS, pH 7.4		

海南省迈迪卡门德商贸有限责任公司

地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502

邮编: 572000

电话: +86-898-31320087

传真: +86-898-31320087

网址: [www.medical-ivd.com](http://www.medical-ivd.com)

邮箱: [info@medical-ivd.com](mailto:info@medical-ivd.com)

Hainan Medikamende Trading Co., Ltd.

Address: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

Phone: +86-898-31320087

Fax: +86-898-31320087

Website: [www.medical-ivd.com](http://www.medical-ivd.com)

Emai: [linfo@medical-ivd.com](mailto:linfo@medical-ivd.com)

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

**Applications:** Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. The following pair is recommended for sandwich immunoassay:

<u>Capture</u>	<u>Detection</u>
L3I00401	L3I00402, L3I00403

**Storage& Shelf life:** Store at  $-20^{\circ}\text{C}$ . Avoid multiple freeze/thaw cycles.

海南省迈迪卡门德商贸有限责任公司

地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502

邮编: 572000

电话: +86-898-31320087 传真: +86-898-31320087

网址: [www.medical-ivd.com](http://www.medical-ivd.com)

邮箱: [info@medical-ivd.com](mailto:info@medical-ivd.com)

Hainan Medikamende Trading Co., Ltd.

Adress: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

Phone: +86-898-31320087 Fax: +86-898-31320087

Website: [www.medical-ivd.com](http://www.medical-ivd.com)

Emai: [linfo@medical-ivd.com](mailto:linfo@medical-ivd.com)