



SPECIFICATION SHEET

Product Name:	Mouse anti human IgG (H&L)McAb	Catalog #:	L1100301
Description:	Mouse Antibody to human Immunoglobulin G (IgG), Heavy & Light chains		
Specificity:	Reacts with heavy chains on human IgG and with the light chains on all human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal human Serum.		
Host Animal:	Mouse		
Immunogen:	Purified human IgG		
Format:	Affinity Purified, Liquid		
Purification:	>90% pure (SDS-PAGE). Immunoaffinity purified. Product is 0.2 micron filtered.		
Concentration:	5.3mg/ml	Lot:	L180601
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY