



SPECIFICATION SHEET

Product Name:	Rabbit anti FITC	Catalog #:	L1101402
Description:	Rabbit Antibody to fluorescein isothiocyanate, FITC PcAb		
Specificity:	Reacts with fluorescein isothiocyanate, FITC		
Host Animal:	Rabbit		
Immunogen:	Purified BSA conjugated FITC		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE). Purified.		
Concentration:	1.2mg/ml	Lot:	L180601
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY