

## SPECIFICATION SHEET

<b>Product Name:</b>	Protein A	<b>Catalog #:</b>	L1I01601
<b>Description:</b>	Recombinant Protein A (Staphylococcus aureus), PA		
<b>Specificity:</b>	Reacts with Fc fragment of IgGs.		
<b>Source:</b>	E.coli,		
<b>Format:</b>	Soluble, recombinant, expressed in E.coli, aqueous solution,		
<b>Purification:</b>	>90% pure (SDS-PAGE). Purified.		
<b>Concentration:</b>	1.74mg/ml	<b>Lot:</b>	C151201
<b>Buffer:</b>	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Use to conjugate with HRP, FITC, Biotin for ELISA, Western blotting and immunohistochemistry, etc.		
<b>Storage&amp; Shelf life:</b>	Store at -20℃ for 3 years. Avoid multiple freeze/thaw cycles.		

**FOR RESEARCH USE ONLY**