



SPECIFICATION SHEET

Product Name: Rabbit anti Chicken IgY **Catalog #:** L1102001

Description: Rabbit anti Chicken Immunoglobulin Y Polyclonal antibody

Specificity: Reacts with heavy chains on Immunoglobulin Y (abbreviated as IgY).

Isotype: IgG

Source: Rabbit

Format: 0.01M PBS with 0.15M NaCl pH 7.4

Concentration: 29mg/ml **Lot:** L190301

Purification: Purified

Application: User optimized

Preservative: None

Storage & Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY