

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Catalog #:</b>	L2G00901	<b>Lot #:</b>	L180901
<b>Description:</b>	Purified Fetal Fibronectin protein		
<b>Background:</b>	FF N, a molecular weight of 310 ~ 550 KD, composed of two polypeptide chains, containing a specific epitope (m CS zone) in which the fetus tumor area.		
<b>Source:</b>	Human placenta		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	About $\geq 90\%$ pure. Affinity purified.		
<b>Concentration:</b>	2 mg/ml		
<b>Buffer:</b>	20mmNH <sub>4</sub> AC PH5.0		
<b>Storage:</b>	Aliquot and store at -20°C. Avoid repeated Freeze/Thaw cycles		

#### Biohazard Information:

Donor samples were tested and found to be negative for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II. Nonetheless caution should be used when handling this material as there is a margin of error in all tests.