



---

SPECIFICATION SHEET

<b>Catalog #:</b>	L4I01102	<b>Lot #:</b>	L180701
<b>Description:</b>	Purified RBP-4 ( Retinol-Binding Protein 4,RBP 4)		
<b>Source:</b>	Human, Native		
<b>Format:</b>	Purified, Lyophilized		
<b>Purification:</b>	>90% (determined by SDS gelelectrophoresis).		
<b>Concentration:</b>	Lypholized,Reconstitute by adding 1ml DD water to 1.0mg/ml		
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Storage:</b>	Store at -20°C. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.