

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Catalog #:</b>	L1V00103	<b>Lot #:</b>	L180301
<b>Description:</b>	Rabbit anti-25-OH-Vitamin D3-BSA (Calcifediol-BSA) Polyclonal Antibody		
<b>Specificity:</b>	Recognizes an epitope on the 25-OH-Vitamin D3-BSA (Calcifediol-BSA)		
<b>Immunogen:</b>	25-OH-Vitamin D3-BSA (Calcifediol-BSA)		
<b>Host Animal:</b>	Rabbit		
<b>Source:</b>	Serum		
<b>Format:</b>	Purified, Liquid		
<b>Concentration:</b>	3.3mg/ml		
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	ELISA		
<b>Storage:</b>	Store at -20°C for long term storage. Avoid repeated freeze-thaw cycles.		