



---

**SPECIFICATION SHEET**

<b>Product Name:</b>	Rabbit anti mouse IgG (H&L)	<b>Catalog #:</b>	L1100901
<b>Description:</b>	Rabbit Antibody to mouse Immunoglobulin G (IgG), Heavy & Light chains		
<b>Specificity:</b>	Reacts with heavy chains on mouse IgG and with the light chains on all mouse immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal mouse Serum.		
<b>Host Animal:</b>	Rabbit		
<b>Immunogen:</b>	Purified Mouse IgG		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure (SDS-PAGE). Purified.		
<b>Concentration:</b>	16.38mg/ml	<b>Lot:</b>	L181101
<b>Buffer:</b>	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	For flow cytometry and fluorescence immunohisto/cytochemistry. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

**FOR RESEARCH USE ONLY**