

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

| | | | |
|-----------------------|--|-------------------|----------|
| ProductName: | Rabbit anti Folic acid (Folate)-BSA | Catalog #: | L6B01103 |
| Synonyms: | Rabbit anti-BSA conjugated FA | | |
| Specificity: | Recognizes an epitope on the Folic acid (Folate)-BSA | | |
| Immunogen: | Folic acid (Folate)-BSA | | |
| Host Animal: | Rabbit | | |
| Source: | Serum | | |
| Format: | Affinity Purified, Liquid | | |
| Concentration: | 5.4mg/ml | Lot: | L181101 |
| Buffer: | 0.01M PBS, pH 7.4 | | |
| Preservative: | None | | |

Description: Folate, distinct forms of which are known as folic acid, folacin, and vitamin B₉, is one of the B vitamins. Synthetic vitamin B₉. Tetrahydrofolate precursor. Promotes DNA synthesis and cell division. Shows anti-inflammatory and chemopreventive effects *in vivo*.

Storage&Shelflife: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.