



SPECIFICATION SHEET

<b>Product Name:</b>	Rabbit anti Hp PcAb	<b>Catalog #:</b>	L1H01202
<b>Description:</b>	Rabbit Antibody to Helicobacter pylori Polyclonal Antibody		
<b>Specificity:</b>	Helicobacter pylori		
<b>Host Animal:</b>	Rabbit		
<b>Immunogen:</b>	Helicobacter pylori		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	>90% pure. Sodium sulfate precipitation and ion-exchange chromatography		
<b>Concentration&amp;Lot:</b>	12mg/ml,	<b>LOT:</b>	L181101
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA, IFA and Western blot. Also suitable for conjugation purposes. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.