

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Description:** Recombinant human MIS/AMH protein,N-terminal .

**Catalog #:** L2F00603

**Description:** Recombinant human anti-Mullerian hormone, AMH

**Source:** E. Coli **Species:** Homo sapiens

**Molecular weight:** 77.7kd (with Tag protein)

**Concentration:** 1.54mg/ml **Lot #:** L181201

**Purity:** >90% by SDS-PAGE

**Formulation:** Liquid

**Application:** ELISA; Other applications have not been tested.  
The optimal dilutions should be determined by end user.

**Storage buffer:** 20mMTrisHCl,PH=8.0, 1mm EDTA,5%Trehalose

**Storage & Expiration:** Upon receipt, aliquot and store at <-20°C  
for long term.Avoid repeated freeze and thaw cycles.

海南省迈迪卡门德商贸有限责任公司  
地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502  
邮编: 572000

Hainan Medikamende Trading Co., Ltd.  
Adress: 6-A-502, Baoyuanhuayuan,  
Xiuyingdadao-No.18, Xiuyingqu,  
Haikou, 572000, Hainan, P.R. China

---

电话: +86-898-31320087 传真: +86-898-31320087 Phone: +86-898-31320087 Fax: +86-898-31320087  
网址: www.medical-ivd.com 邮箱: info@medical-ivd.com Website: www.medical-ivd.com Emai: linfo@medical-ivd.com

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### **Background:**

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.

海南省迈迪卡门德商贸有限责任公司  
地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502  
邮编: 572000

**Hainan Medikamende Trading Co., Ltd.**  
Address: 6-A-502, Baoyuanhuayuan,  
Xiuyingdadao-No.18, Xiuyingqu,  
Haikou, 572000, Hainan, P.R. China

---

电话: +86-898-31320087 传真: +86-898-31320087 Phone: +86-898-31320087 Fax: +86-898-31320087  
网址: www.medical-ivd.com 邮箱: info@medical-ivd.com Website: www.medical-ivd.com Email: linfo@medical-ivd.com