

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Description:** Recombinant human anti-Mullerian hormone Teminal N and C Complex

**Catalog #:** L2F00604

**Description:** Coupling Recombinant human anti-Mullerian hormone Teminal N and C Complex

**Source:** E. Coli      **Species:** Homo sapiens

**Concentration:** 0.2mg/ml (Activation)      **Lot #:** L190101

**Purity:** >90%      by SDS-PAGE

**Formulation:** Liquid

**Application:** ELISA; Other applications have not been tested.  
The optimal dilutions should be determined by end user.

**Storage buffer:** 10mM CBS PH=9.0, 1mm EDTA,5%Trehalose

**Storage & Expiration:** Upon receipt, aliquot and store at <-20°C  
for long term. Avoid repeated freeze and thaw cycles.

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### **Background:**

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.