

SPECIFICATION SHEET

Product Name:	B19 VP1+2	Catalog #:	L4H002
Description:	Recombinant human B19 VP1+VP2 Antigen (B19 VP1+VP2)		
Molecular Weight:	33KD (calculated from sequence) 33KD (determined by SDS gelelectrophoresis)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	> 90% (determined by SDS page)		
Concentration&Lot:	2.5 mg/ml	Lot :	L170701
Buffer:	0.01M PBS pH 7.4 ± 0.2 0.5mM EDTA		
Preservative:	None		
Application:	Immunoassays as EIA, IFA or Immunization.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		