

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	GAD230	Catalog #:	L4I00303
Description:	Recombinant human Glutamic Acid Decarboxylase(GAD230)		
Molecular Weight:	45kd (calculated from sequence) 45KD (determined by SDS gelelectrophoresis)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	> 85% (determined by SDS gelelectrophoresis)		
Concentration& Lot.:	1.47mg/ml (OD280) (Lot. specific)		
Buffer:	Tris-HCL, pH 7.4 ± 0.4		
Preservative:	0.02%NaN3		
Application:	Immunoassays Immunization.		
Storage&Shelflife:			

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.