

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**ProductName:** GAD **Catalog #:** L4I00302

**Description:** Recombinant human Glutamic Acid Decarboxylase(GAD65) Antigen

**Molecular Weight:** 80kd (calculated from sequence)  
80KD(determined by SDS gelelectrophoresis)

**Source:** E. coli

**Format:** Purified, Liquid

**Purification:** > 85% (determined by SDS gelelectrophoresis)

**Concentration& Lot.:** 1.5mg/ml (OD280) **Lot:** L180801

**Buffer:** Tris-HCL, pH 7.4 ± 0.2

**Preservative:** 0.02%NaN3

**Application:** Immunoassays Immunization.

**Storage&Shelflife:**

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.