

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	GAD65	Catalog #:	L4I00304
Description:	Recombinant human Glutamic Acid Decarboxylase(GAD2/GAD65)		
Molecular Weight:	65kd (calculated from sequence) 70KD(determined by SDS gelelectrophoresis)		
Source:	Pichia pastoris		
Format:	Purified, Liquid		
Purification:	> 85% (determined by SDS gelelectrophoresis)		
Concentration& Lot.:	1mg/ml	Lot:	L151201
Buffer:	50mM Tris, 5% glycerol pH:8.0		
Preservative:	None		
Application:	Immunoassays Immunization.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple		