

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName: ICA230 **Catalog #:** L4I00402

Description: Recombinant human Islet cells Antigen(ICA 230)

Molecular Weight: 53kd (calculated from sequence)
53KD(determined by SDS gelelectrophoresis)

Source: E. coli

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelelectrophoresis)

Concentration& Lot.: 0.4mg/ml (Lot. specific)

Buffer: Tris-HCL, pH 7.4 ± 0.2

Preservative: 0.02%NaN3

Application: Immunoassays Immunization.

Storage&Shelflife: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.