

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Description: Recombinant human MIS/AMH protein. **Catalog #:** L2F00602

Description: Recombinant human anti-Mullerian hormone, AMH

Source: E. Coli

Molecular weight: 15 KD

Concentration: 2 mg/ml **Lot #:** L170401

Purity: >90% by SDS-PAGE

Formulation: Liquid

Application: ELISA; Other applications have not been tested.
The optimal dilutions should be determined by end user.

Storage buffer: 50mMTris、PH=8.0, 1mm EDTA,1mMDTT,0.05% SKL

Storage & Expiration: Upon receipt, aliquot and store at <-20°C

for long term.Avoid repeated freeze and thaw cycles.

Background:

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.