

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name: Recombinant Human Placental Growth Factor Protein

Catalog #: L2F01001

SYNONYMS: PGFL; PIGF; PLGF; PIGF-2; D12S1900; SHGC-10760

Source: E. Coli

Species: Homo sapiens

Molecular weight: 19.2KD

Concentration: 1.3mg/ml

Lot #: L20210801

Purity : >90% by SDS-PAGE

Formulation : The protein solution contains 20mM Tris-HCl buffer (pH9.6) ,
5% Trehalose, 1mMEDTA.

Application: ELISA; Positive Control; Immunogen; WB; Activity protein.

Biological Activity: Measured in a cell proliferation assay using
MDA - MB – 231 Human breast cancer cells.

Storage & Expiration: Upon receipt, aliquot and store at <-20°C.

for long term. Avoid repeated freeze and thaw cycles.