

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Description: Recombinant human TSHR protein. **Catalog #:** L2G00601

Description: Thyrotropin receptor, Thyroid-stimulating hormone receptor, TSH-R,

Source: E. Coli

Molecular weight: 80 KD (with Tag protein)

Concentration: 1.6 mg/ml **Lot #:** L170501

Purity: >85% by SDS-PAGE

Formulation: Liquid

Application: ELISA; Other applications have not been tested.
The optimal dilutions should be determined by end user.

Storage buffer: 10mM PBS (pH7.4 ± 0.2)

Storage & Expiration: Upon receipt, aliquot and store at <-20°C
for long term. Avoid repeated freeze and thaw cycles.

海南省迈迪卡门德商贸有限责任公司
地址: 海南省海口市秀英区秀英大道18号宝源花园6-A-502
邮编: 572000
电话: +86-898-31320087, 传真: +86-898-31320087
网址: www.medical-ivd.com 邮箱: info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.
Adress: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,
Xiuyingqu, Haikou, 572000, Hainan, P.R. China
Phone: +86-898-31320087 Fax: +86-898-31320087
Website: www.medical-ivd.com Emai: linfo@medical-ivd.com

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Background:

Thyroid stimulating hormone receptor(TSHR), also known as Thyrotropin receptor is a member of the leucine-richrepeat-containing G protein-coupled receptors. G protein-coupled receptors(GPCRs) are the largest family of plasma membrane receptors. They mediate the effects of several endogenous cues and serve as important pharmacological targets. TSHR is primarily found on the surface of the thyroid epithelial cells. Both TSHR and its ligand TSH have evolved to acquire specificity, minimize cross-reaction to other glycoprotein hormone receptors, and modulate cognate interaction (and thereby thyrotropic activity).