



---

## SPECIFICATION SHEET

<b>Product Name:</b>	TP47	<b>Catalog #:</b>	L2H00702
<b>Description:</b>	<b>Recombinant syphilis TP47 antigen,</b> Antigen Grade		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Protein solution in Tris-Glycin pH8.0		
<b>Purification:</b>	Purity >90%, 47KD		
<b>Concentration&amp;Lot:</b>	1~5mg/ml, lot specific		
<b>Buffer:</b>	20Mm Tris-HCl, pH 8.0 ± 0.02,		
<b>Preservative:</b>	None		
<b>Application:</b>	Immunoassays as EIA, IFA or Immunization.		
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	For research only All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.		