



SPECIFICATION SHEET

Product Name: Recombinant TB-16&38&esat-6 **Catalog #:** L1H01303

Description: Recombinant TB -16,TB-38 and esat-6 antigen,
Antigen Grade

Source: E. coli

Format: Purified, Protein solution in Tris-HCl

Purification: Purity >90%

Concentration&Lot: 2.9 mg/ml **Lot:** L170901

Buffer: 20Mm Tris-HCl, pH 8.0 ± 0.02,

Preservative: None

Applications: EIA,CLIA,IP,ICH,WB

Storage&Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

Warnings: For research only
All cultures have tested negative for HBsAg, HCV, HIV-1 and HIV-2. No test guarantees a product to be non-infectious. Therefore, all material should be considered as potentially infectious.

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.

地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
, Baoyuanhuayuan, Xiuyingdadao-No. 18,
邮编：572000

Address: 6-A-502

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

电话：+86-898-31320087, 传真：+86-898-31320087

Phone: +

86-898-31320087

Fax: +86-898-31320087

网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Website:

www.medical-ivd.com

Email: info@medical-ivd.com