



SPECIFICATION SHEET

Product Name:	TRecombinant TB-16&38	Catalog #:	L1H01302
Description:	Recombinant TB-16 and TB-38 antigen, Antigen Grade		
Source:	E. coli		
Format:	Purified, Protein solution in Tris-HCl		
Purification:	Purity >90%		
Concentration&Lot:	5mg/ml	Lot:	L161201
Buffer:	20Mm Tris-HCl, pH 8.0 ± 0.02,		
Preservative:	None		
Applications:	EIA,CLIA,IP,ICH,WB		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		

海南省迈迪卡门德商贸有限责任公司

Medikamende Trading Co., Ltd.

地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
Baoyuanhuayuan, Xiuyingdadao-No. 18,
邮编：572000

Hainan

Address: 6-A-502,

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

电话：+86-898-31320087, 传真：+86-898-31320087

Phone: +86-898-

31320087 Fax: +86-898-31320087

网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Website: www.medical

-ivd.com Email: info@medical-ivd.com