



---

## SPECIFICATION SHEET

<b>Product Name:</b>	<b>TB CFP-10</b>	<b>Catalog #:</b>	L1H01304
<b>Description:</b>	Recombinant TB CFP-10 antigen, Antigen Grade		
<b>Source:</b>	Recombinant		
<b>Format:</b>	Purified, Protein solution in PB		
<b>Purification:</b>	Purity >90%		
<b>Concentration&amp;Lot:</b>	2.5mg/ml	<b>Lot:</b>	L180501
<b>Buffer:</b>	10Mm PB, pH 7.4 ± 0.02		
<b>Preservative:</b>	None		
<b>Applications:</b>	EIA,CLIA,IP,IHC,WB		
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	For research only		