
SPECIFICATION SHEET

Product Name:	Recombinant TB esat-6	Catalog #:	L1H01302
Description:	Recombinant TB esat-6 antigen, Antigen Grade		
Source:	Recombinant,		
Format:	Purified, Protein solution in Tris-HCl		
Purification:	Purity >90%		
Concentration&Lot:	2.6 mg/ml,	Lot:	L20201201
Buffer:	20Mm Tris-HCl, pH 8.5 ± 0.02,		
Preservative:	None		
Applications:	EIA,CLIA,IP,IHC,WB		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		