

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Description: Recombinant human TSH protein.

Catalog #: L2G00102

Description: Recombinant Thyroid Stimulating Hormone,

Source: E. Coli

Molecular weight: 48 KD (with Tag protein)

Concentration: 1.5mg/ml **Lot #:** L180901

Purity: >85% by SDS-PAGE

Formulation: Liquid

Application: ELISA; Other applications have not been tested.
The optimal dilutions should be determined by end user.

Storage buffer: 10mM PBS (pH7.4 ± 0.2)

Storage & Expiration: Upon receipt, aliquot and store at <-20°C

for long term. Avoid repeated freeze and thaw cycles.