

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Product Name:</b>	Rubella Virus E1 protein, Recombinant (RV-Ag)
<b>Catalog #:</b>	L1Y00101
<b>Description:</b>	Contains a mosaic of the Rubella virus E1 regions. Molecular weight 43kDa. Immunoreactive with Rubella virus positive serum.
<b>Source:</b>	E. coli
<b>Format:</b>	Purified, Liquid
<b>Purification:</b>	>85% pure (SDS PAGE, coomassie staining).
<b>Concentration&amp; Lot.:</b>	4.5mg/ml (Lot: L160701)
<b>Buffer:</b>	10Mm PBS PH7.4 0.5M EDTA
<b>Preservative:</b>	None
<b>Applications:</b>	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.