

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #: L2K002

Lot #: L190101

Description: S-(5'-Adenosyl)-L-homocysteine is S - adenosine methionine dependence methyltransferase reaction in the process of competitive inhibition. It is also known as 5'-Deoxy-S-adenosyl-L-homocysteine AdoHcy ;
S-(5'-Deoxyadenosine-5')-L-homocysteine .

Appearance: White or almost white crystalline powder; odorless; tasteless.

CAS: 979-92-0

Formula: $C_{14}H_{20}N_6O_5S$

MW: 384.41

Form:

Crystalline, Redissolve by 1M HCl:soluble 19.6-20.40mg/ml ,
clear to slightly hazy,colorless to faintly yellow

Contaminants: not detectable

Storage: -20°C or lower, for long term storage.