

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Product Name:** Tox-Ag P30      **Catalog #:** L1Y00301

**Description:** Recombinant T. gondii P30 Ag (P30,SAG1)

**Source:** E.Coli, 30KD  
Strain: RH from ATCC

**Format:** Lysate, Liquid

**Purification:** >85% (SDS-PAGE)

**Concentration&Lot.:** 2mg/ml (Lot. L161001)

**Buffer:** PBS, 0.15M NaCl, pH 7.4 ± 0.2

**Preservative:** None

**Applications:** Specific methodologies have not been tested using this product.

**Storage&Shelf life:**

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.