
SPECIFICATION SHEET

Product name Trypsinogen-2 (TG2) GST-tag

Protein description Recombinant full length Human Trypsinogen-2 (amino acids 1-247) with N terminal GST tag;

Molecular weight 54.5 KD

Protein length 487 aa

Expression host E. coli **Purity** 80% by SDS-PAGE

Purification notes Sepharose derived purification.

Biological activity Immunoreactive with sera of acute pancreatitis, and cystic fibrosis.

Form Liquid

Storage instructions Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Storage buffer

Preservative None

Application ELISA ; WESTERN